

Competitor Fact Sheet: β Beta-bsm[®] & γ Gamma-bsm[®]



For Biocomposites' employees and authorized distributors internal use only.



Feature	genex [®]	β Beta-bsm [®] Injectable γ Gamma-bsm [®] Moldable Putty Zimmer Biomet (Manufactured by Etex [®])
Composition	Powder: 50% β -tricalcium phosphate, 50% CaSO ₄ hemihydrate Mixing solution: Sterile water ¹	Powder: 50% α -tricalcium phosphate, 50% Dicalcium phosphate dihydrate Mixing solution: 0.9% Sodium chloride solution ⁴⁻⁹
Scaffold type	Osteoconductive with negative surface charge for accelerated bone restoration ²	Osteoconductive ^{4,6}
Available sizes	5cc, 10cc	2.5cc, 5cc, 10cc ⁴⁻⁷
Setting time	15 minutes ¹	Isothermally sets hard in 3-5 minutes ^{8,9}
Temperature sensitive setting	No ¹	Yes ^{4,6}
Drillable when fully set	Yes ¹	Yes - during or after setting ^{4,6}
Versatility	Moldable, packable, injectable ¹	<ul style="list-style-type: none"> Beta-bsm: Injectable⁴ Gamma-bsm: Moldable⁶
Injection flexibility	<ul style="list-style-type: none"> Luer Lock syringe with narrow plastic cannula included for hard-to-reach defects OsteoPrecision[™] Graft Delivery Device available to withstand insertion pressure 	<ul style="list-style-type: none"> Beta-bsm: Syringe with 16g needle included⁴ Gamma-bsm: Not injectable
Impurities	No ¹	Unknown

Feature	genex [®]	βBeta-bsm [®] Injectable γGamma-bsm [®] Moldable Putty Zimmer Biomet (Manufactured by Etex [®])
Claimed absorption rate	Up to 12 months ³	Not specified. Expected to be >6 months ¹
Fully absorbs	Yes ¹	May not fully absorb ¹
Dry compressive strength	15MPa ¹	<ul style="list-style-type: none"> • Beta-bsm: 30 MPa⁴ • Gamma-bsm: 35 MPa⁶
Radiopaque	Yes ¹	Yes ¹²
Key selling points and weaknesses	<p>(+) Precisely balanced β-tricalcium phosphate/calcium sulfate hemihydrate¹</p> <p>(+) Contains no Hydroxyapatite (HA) or insoluble impurities¹</p> <p>(+) Fully absorbed within 12 months³</p> <p>(+) Not temperature sensitive¹</p> <p>(+) Provides flexibility with options for molding, packing and injecting</p> <p>(+) Drillable when fully set¹</p> <p>(+) Radiopaque¹</p> <p>(+) Negatively charged surface chemistry accelerates bone growth up to 5x normal levels²</p> <p>(+) Restores bone to normal trabecular structure in 36 weeks³</p> <p>(+) Comprehensive support network for our customers and hospitals</p>	<p>(+) Gamma-bsm putty hardens in a wet environment^{4,6}</p> <p>(+) Can be drilled during or after setting^{4,6}</p> <p>(-) Calcium phosphate converts to a poorly crystalline HA during setting^{1,4}</p> <p>(-) HA has a slow and incomplete absorption rate¹</p> <p>(-) HA can cause a long-term nidus for infection¹</p> <p>(-) Contraindicated for infected sites^{10,11}</p> <p>(-) Temperature sensitive - setting time is an inverse function of local body temperature^{4,6,7,10,11}</p> <p>(-) Must use different products when injecting or molding</p> <p>* (+) = competitor selling points (-) = competitor weaknesses</p>

References:

1. Biocomposites, Data on file.
2. Cooper, J.J., J.A. Hunt, and F. Pu, Enhancing the Osteogenic Potential of Bioabsorbable Implants through Control of Surface Charge. Presented at the Society for Biomaterials 2007 Annual Meeting. 2007: Chicago, Illinois, USA.
3. Yang HL et al. Bone healing response to a synthetic calcium sulfate/beta-tricalcium phosphate graft material in a sheep vertebral body defect model. J Biomed Mater Res B Appl Biomater. 2012;100B(7):1911-21.
4. Zimmer Biomet βBeta-bsm[®] Injectable Bone Substitute Material Flyer. 1671.1-US-en. 2017.
5. Zimmer Biomet βBeta-bsm[®] Injectable Bone Substitute Material Mixing Chart. 0919.1-GLBL-en-REVO718. 2018.
6. Zimmer Biomet γGamma-bsm[®] Moldable Putty Bone Substitute Material Flyer. 1672.1-US-en. 2017.
7. Zimmer Biomet γGamma-bsm[®] Moldable Putty Bone Substitute Material Mixing Chart. 0918.1-GLBL-en-REV1018. 2018.
8. <https://www.zimmerbiomet.com/medical-professionals/biologics/product/beta-bsm-injectable-bone-substitute-material.html>
9. <https://www.zimmerbiomet.com/medical-professionals/biologics/product/gamma-bsm-moldable-putty.html>
10. ETEx Corporation, Zimmer GmbH βBeta-bsm[®] Injectable Bone Substitute Material Instructions for Use. 50-1478v05.
11. ETEx Corporation, Zimmer GmbH γGamma-bsm[™] Moldable Bone Substitute Material Instructions for Use. 50-144lv04.
12. Russell T.A., Leighton R.K. and on behalf of the Alpha-BSM Tibial Plateau Fracture Study Group, Comparison of Autogenous Bone Graft and Endothermic Calcium Phosphate Cement for Defect Augmentation in Tibial Plateau Fractures. A Multicenter, Prospective, Randomized Study. J Bone Joint Surg Am. 2008;90:2057-2061.

For indications, contraindications, warnings and precautions see Instructions for Use.

©2020, Biocomposites, genex and Power to Restore are trademarks/registered trademarks of Biocomposites Ltd. All rights reserved.
No unauthorized copying, reproduction, distributing or re-publication is allowed unless prior written permission is granted by the owner, Biocomposites Ltd.

Patents granted: EP I390086 BI, US 8632796, CN ZL02809194.9, US 8496955

MA0332R1