

Competitor Fact Sheet: HydroSet® & HydroSet® XT



For Biocomposites' employees and authorized distributors internal use only.



Feature	genex®	HydroSet® & HydroSet® XT Stryker®
Composition	Powder: 50% β -tricalcium phosphate, 50% CaSO_4 hemihydrate Mixing solution: Sterile water ¹	Powder: Dicalcium phosphate dihydrate, Tetracalcium phosphate, Tri-sodium citrate Mixing solution: Sodium phosphate, Polyvinylpyrrolidone, water ⁴
Scaffold type	Osteoconductive with negative surface charge for accelerated bone restoration ²	Osteoconductive ^{4,5,6}
Available sizes	5cc, 10cc	3cc, 5cc, 10cc, 15cc ^{6,7}
Setting time	15 minutes ¹	Within 8 minutes when mixed according to manufacturer's instructions ⁷
Temperature sensitive setting	No ¹	Yes ^{4,5,6}
Drillable when fully set	Yes ¹	Yes ^{4,5}
Versatility	Moldable, packable, injectable ¹	Moldable, injectable ⁶
Injection flexibility	<ul style="list-style-type: none"> Luer Lock syringe with narrow plastic cannula included for hard-to-reach defects OsteoPrecision™ Graft Delivery Device available to withstand insertion pressure 	<ul style="list-style-type: none"> Hydroset: Standard delivery syringe with cannula included⁴ Hydroset XT: Pre-filled, self-contained mixing and delivery system with locking torque handle and cannula^{7,8}
Impurities	No ¹	Unknown

Feature	genex [®]	HydroSet [®] & HydroSet [®] XT Stryker
Claimed absorption rate	Up to 12 months ³	Not specified. Expected to be >6 months ¹
Fully absorbs	Yes ¹	No ¹
Dry compressive strength	15 MPa ¹	15.9MPa ⁹
Radiopaque	Yes ¹	Yes ^{6,9}
Key selling points and weaknesses	(+) Precisely balanced β -tricalcium phosphate/calcium sulfate hemihydrate ¹ (+) Contains no HA or insoluble impurities ¹ (+) Fully absorbed within 12 months ³ (+) No contraindication against use in articulating surfaces ¹ (+) Not temperature sensitive ¹ (+) Provides options for injection flexibility (+) Drillable when fully set ¹ (+) Radiopaque ¹ (+) Negatively charged surface chemistry accelerates bone growth up to 5x normal levels ² (+) Restores bone to normal trabecular structure in 36 weeks ³ (+) Comprehensive support network for our customers and hospitals	(+) Hydroset XT features a self-contained delivery system for increased working time ⁸ (+) Drillable ^{4,7} (+) Hardens in a wet environment ^{4,6} (+) Radiopaque ⁶ (-) Calcium phosphate converts to Hydroxyapatite (HA) after implantation ^{4,7} (-) HA has a slow and incomplete absorption rate ¹ (-) HA can cause a long-term nidus for infection ¹ (-) Contraindicated for infected sites ⁵ (-) Contraindicated for bone voids that link joint spaces and articulating surfaces ⁵ (-) Temperature sensitive ⁴ (-) Contact and heat transfer between hands and syringe may decrease injection time ^{4,6} (-) Setting time shown to extend up to 19 minutes ¹⁰ * (+) = competitor selling points (-) = competitor weaknesses

References:

1. Biocomposites, Data on file.
2. Cooper, J.J., J.A. Hunt, and F. Pu, Enhancing the Osteogenic Potential of Bioabsorbable Implants through Control of Surface Charge. Presented at the Society for Biomaterials 2007 Annual Meeting. 2007: Chicago, Illinois, USA.
3. Yang HL et al. Bone healing response to a synthetic calcium sulfate/beta-tricalcium phosphate graft material in a sheep vertebral body defect model. J Biomed Mater Res B Appl Biomater 2012;100B(7):1911-21.
4. Stryker HydroSet™ Injectible HA Bone Substitute brochure. LHS-B MS/GS 3C 08/06. 2006.
5. Stryker HydroSet Injectible HA Bone Substitute Osteosynthesis brochure. 90-07900 Lot B1008. 2008.
6. Stryker Craniomaxillofacial Hydroset Fast-setting HA bone substitute brochure. CMF-BR-148. Rev None 17430. 2018.
7. Stryker Bone graft substitute portfolio. BIAL BR6 Rev 117562. 2018.
8. Stryker Spine HydroSet® XT Injectible Bone Substitute brochure. BIHYD-SS-1 Rev 2 15770. 2017.
9. Stryker HydroSet™ Injectible HA Bone Substitute At-A-Glance Reference Guide. LHSAG-B MS/GS 1M 09/06. 2006.
10. Niall Kent et al. In vitro and in vivo study of commercial calcium phosphate cement HydroSet™. 2016 Wiley Periodicals, Inc. J Biomed Mater Res Part B: Appl Biomater, 106(B):21-30, 2018.

For indications, contraindications, warnings and precautions see Instructions for Use.

©2020, Biocomposites, genex and Power to Restore are trademarks/registered trademarks of Biocomposites Ltd. All rights reserved.
No unauthorized copying, reproduction, distributing or re-publication is allowed unless prior written permission is granted by the owner, Biocomposites Ltd.

Patents granted: EP 1390086 B1, US 8632796, CN ZL02809194.9, US 8496955

MA0312RI