

## Versatility at your fingertips

genex is versatile and easy to prepare. It comes as a paste or putty in a range of volumes.

PRODUCT	PASTE VOLUME	SETTING TIMES	ORDER CODES
genex Setting paste for moulding, injecting and packing	5cc	15 minutes approx.	900-005
	10cc		900-010
genex ds Setting paste for difficult-to-reach surgical sites	2.5cc	15 minutes approx.	980-002
	5cc		980-005
genex putty Non-setting putty that requires no mixing	2.5cc	Non-setting	920-002
	5cc		920-005
	10cc		920-010

Personal note

For indications, contraindications, warnings and precautions see Instructions for Use.

©2018, Biocomposites, genex and Power to Restore are trademarks / registered trademarks of Biocomposites Ltd. All rights reserved. No unauthorised copying, reproduction, distributing or republication is allowed unless prior written permission is granted by the owner, Biocomposites Ltd.

Patents granted: EP I390086 BI, US 8632796, CN ZL02809194.9, US 8496955

**References:** **1.** Yang, H.L., et al, Bone healing response to a synthetic calcium sulfate / beta-tricalcium phosphate graft material in a sheep vertebral body defect model. J Biomed Mater Res B Appl Biomater, 2012. 100B(7): p. 1911-21. **2.** Data on File, Mr A Nisar and Mr S Gopal. **3.** Data on File, Mr HK Sharma. **4.** Data on File, Prof JB Richardson.

POWER TO RESTORE WITHOUT LEAVING A TRACE



The only remaining  
evidence of the trauma

## Perfect partner for your trauma and non-unions

geneX is a catalyst for bone healing. It complements the body's natural healing processes and encourages normal bone structure to be restored at a steady rate. Over 12 months, it is completely absorbed and leaves no foreign artefacts to impair structural integrity.<sup>1,2</sup>

- ✓ Approved for bony voids and defects that are not intrinsic to structural stability

geneX offers successful long-term outcomes across a range of surgical applications, including tibial plateau fractures and long-bone non-unions<sup>2-4</sup>



TIBIAL PLATEAU FRACTURES<sup>3</sup>



PROXIMAL FEMUR FRACTURES<sup>2</sup>



DISTAL TIBIA NON-UNIONS<sup>4</sup>

## Designed to be completely absorbed and leave no trace

geneX is a precisely balanced  $\beta$ -tricalcium phosphate/calcium sulfate hemihydrate compound with distinct design properties.

- ✓ Contains no hydroxyapatite (HA)
- ✓ Fully absorbed within 12 months
- ✓ Enhances osteogenic response
- ✓ Restores strong healthy bone



Discover the  
true potential at  
**biocomposites.com**