

STIMULAN®

Case-by-case flexibility

Can be used in three bead sizes, as a paste or for injection,
STIMULAN optimizes packing in bone voids and defects.

PRODUCT	PASTE VOLUME	BEAD VOLUME	SETTING TIMES	ORDER CODES
STIMULAN Rapid Cure Rapid setting paste with individual bead kit	5cc	12cc	Approx. 4 minutes	620-005
	10cc	25cc		620-010
	20cc	50cc		620-020
STIMULAN KIT Standard setting, fully injectable paste with individual bead kit	5cc	10cc	Approx. 8 minutes	600-005
	10cc	20cc		600-010

Personal note

References: 1. Biocomposites Ltd, STIMULAN, Instructions for Use (US). 2. Data on file, Mr. Sharma 3. Data on file, Dr. Schlatterer 4. Data on file, Dr. Webb 5. Data on file.

For indications, contraindications, warnings and precautions see Instructions for Use. Concurrent use of locally administered antibiotics may affect setting time.

This brochure may include the use of STIMULAN or techniques that go beyond the current clearance / approval granted by the relevant regulatory authority. Please contact your local representative for further information.

©2018, Biocomposites, STIMULAN, Bringing Calcium to Life and Power to Transform are trademarks/registered trademarks of Biocomposites Ltd.

All rights reserved. No unauthorized copying, reproduction, distributing or republication is allowed unless prior written permission is granted by the owner, Biocomposites Ltd.

Patents granted: GB2367552, EP 1204599 B1, US 6780391, EP 2594231 B1, US 8883063, CN ZL201210466117.X, GB2496710, EP 3058899 B1

Patents pending: GB1502655.2, US 15/040075, CN 201610089710.5, US 15/288328, GB1704688.9



STIMULAN[®]

Perfect partner for cases with an infected site

STIMULAN is a truly absorbable calcium sulfate, specifically designed to complement your dead space management, across a wide range of settings.¹

- ✓ Cleared for placement directly at site of infection in bone voids and defects
- ✓ May help improve outcomes and reduce associated costs of care


STIMULAN can be used in cases with an infected site to fill a bone void or defect created by a cyst, tumor, traumatic injury or osteomyelitis¹



NON-UNION²



OSTEOMYELITIS³



TRAUMA⁴

Recrystallized for consistent and reliable performance

Only STIMULAN undergoes a patented recrystallization process that results in its consistent and reliable performance.^{1,5}

- ✓ Controlled purity
- ✓ Physiologic pH level
- ✓ Truly absorbable at an optimum rate
- ✓ No third body damage



Discover the
true potential at
biocomposites.com